

Rayat Shikshan Sansthas  
**Rajarshi Chhatrapati Shahu College Kolhapur**

**Department of Botany**

**WALL PAPER**

<b>Sr No</b>	<b>Topic of Wall paper</b>	<b>Name of Students</b>	<b>Class</b>
1	Detect Adulteration in Milk and Milk Product	Miss. Ware Trupti Krushnat	B.Sc-II
2	Detect Adulteration in Beverages	Miss. Thakur Kajal Bajrang	B.Sc-III
3	Detect Adulteration in Fruits & Vegetables	Mr. Kamble Krishna Laxman	B.Sc-III
4	Detect Adulteration in salts, Spices & Condiments	Mr.Shivputra Kempayya Hanchinwar	B.Sc-III
5	Detect Adulteration in Food Grain & Its Product	Mr. Kamble Pranam Kerba	B.Sc-III
6	Detect Adulteration in Sugar & Confectionery	Mr. Yesare Aniruddha Baban	B.Sc-III
7	Detect Adulteration in Oils & Fats	Miss. Suryawanshi Supriya Tukaram	B.Sc-III

**Detect Adulteration**  
**in**  
**Milk and Milk Product**

Name:Ware Trupati

Class:B.Sc-II

---

## Detection of water in milk

### Testing method:

- 1 Put a drop of milk on a polished slanting surface.
- 2 Pure milk either stays or flows slowly leaving a white trail behind.
- 3 Milk adulterated with water will flow immediately without leaving a mark.



Pure milk



Adulterated milk

---

## Detection of detergent in milk

### Testing method:

- 1 Take 5 to 10ml of sample with an equal amount of water.
- 2 Shake the contents thoroughly.
- 3 If milk is adulterated with detergent, it forms dense lather.
- 4 Pure milk will form very thin foam layer due to agitation..



Pure milk



Adulterated milk

---

## Detection of starch in milk and milk products (khoya, chenna, paneer)

### Testing method:

- 1 Boil 2-3 ml of sample with 5ml of water.
- 2 Cool and add 2-3 drops of tincture of iodine.
- 3 Formation of blue colour indicates the presence of starch.  
(In the case of milk, addition of water and boiling is not required)



Pure milk



Adulterated milk

---

## Detection of mashed potatoes, sweet potatoes and other starches in ghee/butter

### Testing method:

- 1 Take ½ teaspoon of ghee/butter in a transparent glass bowl.
- 2 Add 2-3 drops of tincture of iodine.
- 3 Formation of blue colour indicates the presence of mashed potatoes, sweet potatoes and other starches.



Pure



Adulterated

# **Detect Adulteration in Beverages**

**Name: Thakur Kalal Bajrang**

**Class:B.Sc-III (Botany)**




---

## Detection of clay in coffee powder

### Testing method:

- 1 Add ½ teaspoon of coffee powder in a transparent glass of water.
- 2 Stir for a minute and keep it aside for 5 minutes. Observe the glass at the bottom.
- 3 Pure coffee powder will not leave any clay particles at the bottom.
- 4 If coffee powder is adulterated, clay particles will settle at the bottom.



Pure



Adulterated

---

## Detection of chicory powder in coffee powder

### Testing method:

- 1 Take a transparent glass of water.
- 2 Add a teaspoon of coffee powder.
- 3 Coffee powder floats over the water but chicory begins to sink.



Pure



Adulterated

# Detection of exhausted tea in tea leaves

## Testing method - 1:

- 1 Take a filter paper and spread few tea leaves.
- 2 Sprinkle with water to wet the filter paper.
- 3 Wash the filter paper under tap water and observe the stains against light.
- 4 Pure tea leaves will not stain the filter paper.
- 5 If coal tar is present, it will immediately stain the filter paper.

## Testing method - 2:

- 1 Take small amount of tea leaves/ dust and place it in the centre of a filter paper.
- 2 Add water drop by drop at the heap of the tea leaves/ dust.
- 3 If the tea is adulterated with a coloured tea, water will dissolve the added colour and leave streak of colour in the filter paper.



Pure



Adulterated

## Testing method - 3:

- 1 Spread a little slaked lime on white porcelain tile or glass plate
- 2 Sprinkle a little tea dust on the lime.
- 3 Red, orange or other shades of colour spreading on the lime will show the presence of coal tar colour.
- 4 In case of genuine tea, there will be only a slight greenish yellow colour due to chlorophyll, which will appear after some time.

---

## Detection of iron filings in tea leaves

### Testing method:

- 1 Take small quantity of tea leaves in a glass plate.
- 2 Move the magnet through the tea leaves.
- 3 Pure tea leaves will not show any iron filings on the magnet.
- 4 If adulterated, then iron filings will be seen on the magnet.



Pure



Adulterated



# **Detect Adulteration in Fruits & Vegetables**

**Name:** Kamble Krishna Laxman

**Class:** B.Sc-III (Botany)

## Detection of malachite green in green vegetables like bitter gourd, green chilli and others.

### Testing method - 1:

- 1 Take a cotton piece soaked in water or vegetable oil. (conduct the test separately)
- 2 Rub the outer green surface of a small part of green vegetable/chilli.
- 3 If the cotton turns green, then it is adulterated with malachite green.



Pure



Adulterated

### Testing method - 2:

- 1 Take a small part of the sample and place on a piece of moistened white blotting paper.
- 2 The impression of colour on the paper indicates the use of malachite green, or any other low priced artificial colour.



Green Vegetables



Malachite Green

---

## Detection of artificial colour on green peas

### Testing method:

- 1 Take little amount of green peas in a transparent glass.
- 2 Add water to it and mix well.
- 3 Let it stand for half an hour.
- 4 Clear separation of colour in water indicates adulteration.



Pure



Adulterated

---

## Detection of rhodamine B in sweet potato

### Testing method:

- 1 Take a cotton ball soaked in water or vegetable oil. (conduct the test separately)
- 2 Rub the outer red surface of the sweet potato.
- 3 If cotton absorbs colour, then it indicates the usage of rhodamine B for colouring the outer surface of sweet potato.



Pure



Adulterated

# **Detect Adulteration in Salts, Spices & Condiments**

**Name:** Shivputra Kempayya Hanchinwar

**Class:** B.Sc-III (Botany)



## Detection of foreign resin in asafoetida (hing)

### Testing method-1:

- 1 Burn small quantity of asafoetida in a stainless steel spoon.
- 2 Pure asafoetida will burn like camphor.
- 3 Adulterated asafoetida will not produce bright flame like camphor.



Pure



Adulterated

### Testing method-2:

- 1 Powder a gram of asafoetida and take it in a glass container.
- 2 Add one tea spoon of water. Mix thoroughly by shaking.
- 3 Milky white solution with no sediments represents pure asafoetida.



Asafetida



Non-Edible Gum/Resin

---

## Detection of papaya seeds in black pepper

### Testing method - 1:

- 1 Add some amount of black pepper to a glass of water.
- 2 Pure black pepper settles at the bottom.
- 3 In the adulterated black pepper, papaya seeds float on the surface of water.



**Black pepper**




**Papaya seeds**

## Detection of papaya seeds in black pepper

### Testing method - 2:

- 1 Spread spice on a white paper.
- 2 Observe the appearance of the sample using the magnifying glass.
- 3 Black pepper is brown in colour. It has a wrinkled surface and has a characteristic smell and pungent taste.
- 4 The papaya seeds have shrunken smooth surface and oval shape. It is greenish brown or blackish brown in colour and has a repulsive flavour.



**Black pepper**



**Papaya seeds**

---

## Detection of light black berries in black pepper

### Testing method

- 1 Press the berries with the help of fingers.
- 2 Light berries will break easily while black berries of pepper will not break.



**Black pepper**



**Light black berries**

---

## Detection of soap stone or other earthy matter in asafoetida (hing)

### Testing method:

- 1 Shake little portion of the sample with water and allow to settle.
- 2 Pure asafoetida will not leave any soap stone or other earthy matter at the bottom.
- 3 If asafoetida is adulterated, soap stone or other earthy matter will settle down at the bottom.



**Pure**



**Adulterated**



---

## Detection of artificial/water soluble synthetic colours in chilli powder

### Testing method:

- 1 Sprinkle chilli powder on the surface of water taken in a glass tumbler.
- 2 The artificial colourants will immediately start descending in colour streaks.



Pure



Adulterated

---

## Detection of light black berries in black pepper

### Testing method:

- 1 Float the sample of black pepper in alcohol (rectified spirit).
- 2 The mature black pepper berries sink while the light black pepper floats.



Black pepper



Light black berries



---

## Detection of saw dust in chilli powder

### Testing method:

- 1 Add the sample to water.
- 2 The saw dust will float at the surface of water while Chilli powder will settle down in bottom.



Chilli powder



Saw dust

---

## Detection of starch in asafoetida

### Testing method:

- 1 Tincture of iodine is added to the sample of asafoetida.
- 2 Appearance of blue colour shows the presence of starch.



Asafoetida



Starch

---

## Detection of chalk in common salt

### Testing method:

- 1 Stir a spoonful of sample of salt in a glass of water.
- 2 The presence of chalk will make solution white and other insoluble impurities will settle down.



Common salt



Chalk

---

## Detection of exhausted cloves in cloves

### Testing method:

- 1 Take some water in a glass and put cloves.
- 2 Genuine cloves will settle down at the bottom while exhausted cloves will float on surface.



Cloves



Volatile oil extracted cloves

---

## Detection of cassia bark in cinnamon

### Testing method:

- 1 Take small quantity of cinnamon in a glass plate.
- 2 If adulterated, on close visual examination, cassia bark that comprises of several layers in between the rough outer and inner most smooth layers can be differentiated from cinnamon.
- 3 Cinnamon barks are very thin that can be rolled around a pencil or pen. It has a distinct smell.



Cinnamon



Cassia

---

## Detection of grass seeds coloured with charcoal dust in cumin seeds

### Testing method:

- 1 Rub small amount of cumin seeds on palms.
- 2 If palms turn black, adulteration is indicated.



Pure



Adulterated



---

## Detection of argemone seeds in mustard seeds

### Testing method:

- 1 Take small quantity of mustard seeds in a glass plate.
- 2 Examine visually for the argemone seeds.
- 3 Mustard seeds have a smooth surface and when pressed, inside it is yellow in colour.
- 4 Argemone seeds have grainy, rough surface and are black in colour. When pressed, it is white in colour from inside.



**Mustard seeds**



**Argemone seeds**

---

## Detection of lead chromate in turmeric whole

### Testing method:

- 1 Add small quantity of turmeric whole in a transparent glass of water.
- 2 Pure turmeric will not leave any colour.
- 3 Adulterated turmeric appears to be bright in colour and leaves colour immediately in water.



**Pure**



**Adulterated**



---

## Detection of artificial colour in turmeric powder

### Testing method:

- 1 Add a teaspoon of turmeric powder in a glass of water.
- 2 Natural turmeric powder leaves light yellow colour while settling down.
- 3 Adulterated turmeric powder will leave a strong yellow colour in water while settling down.



Pure



Adulterated

---

## Detection of sawdust and powdered bran in powdered spices

### Testing method:

- 1 Sprinkle powdered spices on the water surface.
- 2 Pure spices will not leave any saw dust/powdered bran on the surface of water.
- 3 If spices are adulterated, saw dust/powdered bran will float on the surface.



Pure



Adulterated

---

## Differentiation of common salt and iodised salt

### Testing method:

- 1 Cut a piece of potato, add salt and wait for a minute.
- 2 Add two drops of lemon juice.
- 3 If it is iodised salt, blue colour will develop.
- 4 In the case of common salt, there will be no blue colour.



**Iodised salt**



**Common salt**

---

## Detection of coloured dried tendrils of maize cob in saffron

### Testing method:

- 1 Genuine saffron will not break easily like artificial. Artificial saffron is prepared by soaking maize cob in sugar and colouring it with coal tar.
- 2 Take a transparent glass of water and add small quantity of saffron.
- 3 If saffron is adulterated, the artificial colour dissolves in water rapidly. A bit of pure saffron when allowed to dissolve in water will continue to give its saffron colour so long as it lasts.



**Saffron**



**Coloured tendrils**

# **Detect Adulteration in Food Grains & Its Products**

Name: Kamble Pranam Kerba

Class:B.Sc-III (Botany)



---

## Detection of extraneous matter (dust, pebble, stone, straw, weed seeds, damaged grain, weeviled grain, insects, rodent hair and excreta) in food grains

### Testing method:

- 1 Take small quantity of sample in a glass plate.
- 2 Examine the impurities visually.
- 3 Pure food grains will not have any such impurities.
- 4 Impurities are observed visually in adulterated food grains.



Pure



Adulterated

---

## Detection of dhatura in food grains

### Testing method:

- 1 Take small quantity of food grains in a glass plate.
- 2 Examine the impurities visually.
- 3 Dhatura seeds which are flat with edges and blackish brown in colour can be separated out by close examination.



Dhatura seeds in food grains



Dhatura seeds



---

## Detection of excess bran in wheat flour

### Testing method:

- 1 Take a transparent glass of water.
- 2 Sprinkle a spoon of wheat flour on the surface of water.
- 3 Pure wheat flour will not show excess bran on water surface.
- 4 Adulterated wheat flour shows excess bran floating on water surface.



Pure wheat flour



Excess bran in wheat flour

---


## Detection of khesari dal in dal whole and split

### Testing method:

- 1 Take small quantity of dal whole or split in a glass plate.
- 2 Examine the impurities visually.
- 3 Khesari dal which has edged type appearance showing a slant on one side and square in appearance can be separated out by close examination.
- 4 Pure dal will not have any such impurities.



Pure dal



Khesari dal

---

## Detection of added colour in food grains

### Testing method:

- 1 Take a transparent glass of water.
- 2 Add 2 teaspoons of food grains and mix thoroughly.
- 3 Pure food grains will not leave any colour.
- 4 Adulterated food grains leaves colour immediately in water.



Pure



Adulterated

---

## Detection of turmeric in sella rice

### Testing method:

- 1 Take a tea spoon of rice in a glass plate.
- 2 Sprinkle a small amount of soaked lime (commonly known as chuna which is used in pan) on the rice grains.
- 3 Pure grains will not form red colour.
- 4 Adulterated grains will form red colour.



Pure



Adulterated



---

## Detection of rhodamine B in ragi

### Testing method:

- 1 Take cotton ball soaked in water or vegetable oil. (conduct the test separately)
- 2 Rub the outer surface of the ragi.
- 3 If cotton absorbs colour, then it indicates the adulteration of rhodamine B for colouring the outer surface of ragi.



Pure



Adulterated

---

## Detection of chakunda beans in pulses

### Testing method:

- 1 Take small quantity of pulses in a transparent glass plate.
- 2 Examine the impurities visually.
- 3 Chakunda beans can be separated out by close examination.



Chakunda beans

---

## Detection of sand, soil, insects, webs, lumps, rodent hair and excreta in Atta, Maida, Suji (Rawa)”

### Testing method:

- 1 These can be identified by visual examination.



**Atta, Maida, Suji (Rawa)**



**Sand, Soil, insects, webs, lumps,  
rodent hair and excreta**



# **Detect Adulteration in Sugars & Confectionery**

**Name: Yesare Aniruddha Baban**

**Class:B.Sc-III (Botany)**

---

## Detection of sugar solution in honey

### Testing method - 1:

- 1 Take a transparent glass of water.
- 2 Add a drop of honey to the glass.
- 3 Pure honey will not disperse in water.
- 4 If the drop of honey disperses in water, it indicates the presence of added sugar.



Pure

### Testing method - 2:

- 1 Take a cotton wick dipped in a pure honey and light with a match stick.
- 2 Pure honey will burn.
- 3 If adulterated, the presence of water will not allow the honey to burn if it does; it will produce a cracking sound.



Adulterated

---

## Detection of chalk powder in sugar/pithi sugar/jaggery

### Testing method:

- 1 Take a transparent glass of water.
- 2 Dissolve 10g of sample in water.
- 3 If sugar/pithi sugar/jaggery is mixed with chalk, the adulterant will settle down at the bottom.



Pure



Adulterated

## Detection of aluminium leaves in silver leaves

### Testing method:

- 1 Take some portion of the leaf and crush it between two fingers.
- 2 Pure silver leaves will be easily crushed and crumble to the powder form while aluminium leaves will only break into smaller shreds.
- 3 Further take the suspected silver leaves and make it in the form of a ball and burn it with the help of a flame.
- 4 Pure silver leaves burn away completely leaving glistening balls while aluminium leaves are reduced to grey ash.



**Silver leaves**



**Aluminium leaves**



**Silver leaves**



**Aluminium leaves**

# **Detect Adulteration in Oils & Fats**

**Name: Suryawanshi Supriya Tukaram**

**Class:B.Sc-III (Botany)**



---

## Detection of other oils in coconut oil

### Testing method:

- 1 Take coconut oil in a transparent glass.
- 2 Place this glass in refrigerator for 30 minutes. (Do not keep in the freezer)
- 3 After refrigeration, coconut oil solidifies.
- 4 If coconut oil is adulterated, then other oils remain as a separate layer.



Pure



Adulterated

---

## Detection of TOCP (Tri-Ortho-Cresyl-Phosphate) in oils and fats

### Testing method:

- 1 Take 2ml of sample of oil.
- 2 Add on a little amount of yellow butter (Solid).
- 3 Immediate formation of red colour indicates the presence of TOCP.



Pure



Adulterated

---

## Proper winterization of refined winterized salad oils

### Testing method:

- 1 Take 100ml sample oil in bottle, cork tightly and seal with paraffin.
- 2 The bottle is completely submerged in bucket containing finely cracked ice and water is added until it rises to top of the bottle.
- 3 The bucket is kept filled solidly with ice by removing any excess water and adding ice when necessary.
- 4 After 5.5hours remove the bottle and examine oil.
- 5 If it is properly winterised, sample will be brilliant, clear and limpid.



**Refined winterized salad oils**